NEBUALDO

Easy Pro

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Mesh Nebulizer

Manual de usuario / User manual

Aldo-Unión
www.aldo-union.com

NEBUALDO

Easy Pro

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Mesh Nebulizer

Manual de usuario / User manual

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Be sure to read this user manual before using the unit in order to use it safely and correctly

1. Device description

Product Instructions

Thanks to self-excited oscillation circuit, the ceramic sheets contained in nebulizing head together with the metal mesh are driven successively to oscillate at a high frequency, which makes the liquid flow through the micromesh of the metal mesh and be ejected and nebulized to form inhalant aerosol. The medicine passes through the inhalation mask to the patient's respiratory system, in order to achieve the purpose of inhalation treatment.


Symbols

The warning signs and the sample icons shown here are listed for you to use this product safely and correctly as well as to prevent the risk and the damage to you and others from happening.

 Manufacturer	 Operating instructions
 The device should be sent to the special agencies according to local regulations for separate collection after its shelf-life	 Type BF Applied Part
 Medical device	 This symbol shall be accompanied by the manufacturer's serial number

Components/Accessories

Before installing and using the product, the user first checks whether the appearance of the product is good, whether the variety and quantity of the accessories are consistent with the list of accessories described in the instructions, and if there is any defect, the user shall contact the supplier or manufacturer immediately.

1. Medicine cup
2. Mask
3. String
4. Main unit
5. Main unit power button
6. Main unit connection cable
7. Type C adapter head

2. Specifications

Product Name	Nebualdo Easy Pro
Model	VP-M10
Power source	DC 5V 1A
Nebulization Rate	≥0.20ml/min
Particle Size	MMAD approximately 5µm
Sound Level	≤50dB
Ultrasonic Frequency	Approximately 114kHz
Capacity of medicine cup	8ml
Water Temperature of medicine cup	≤60°C
Weight of medicine cup	Approximately 11g
Dimensions	42 (W) mm x 34 (H) mm x 40 (D) mm
Safety Classification	Equipment that cannot be used in the case of flammable narcotic gas mixed with air or flammable narcotic gas mixed with oxygen or nitrous oxide
Scope of Application	Atomized drugs for treatment
Contraindications	None
Electrical safety requirements	Internal power supply equipment, BF application part
Operation mode	Continuous operation
Working Environment	Temperature: 10°C to 40°C Humidity: 30 to 85%RH Atmospheric Pressure: 86kPa-106kPa
Transport and Storage Environment	Temperature: -20 to 50°C Humidity: 30 to 85%RH Atmospheric Pressure: 86kPa-106kPa

3. Safety precaution

Warnings:

- ⚠ For the type, dose and regimen of the medication, be sure to follow the instructions of a doctor.
Your symptom(s) may worsen.
- ⚠ Clean and disinfect medicine cup and inhalation mask, after each use.
- ⚠ Be sure to dry the cleaned and disinfected parts promptly, and store them in a clean place.
Miscellaneous bacteria may propagate and you may get infected.
- ⚠ This unit is designed for personal use. It is not appropriate for more than one user.
- ⚠ Keep the unit out of reach of children or infants.
They may swallow small pieces. The main unit may be damaged and in case it is damaged it can cause injury to children or infants.
In case a child swallows small pieces, immediately consult a doctor immediately.
- ⚠ Do not wash the main unit with water or hot water.
Do not immerse the unit in water since it may experience electricity leakage or you may suffer electric shock.
- ⚠ This device is a professional atomization treatment equipment. Children and adults who are unable to take care of themselves should ensure that a guardian is present when using it.

Caution:

- ⚠ Do not inhale by using water in the medicine cup. Your symptom(s) may worsen.
- ⚠ Do not drop medication on the main unit. If you drop medication, wipe it off immediately with gauze. If you use the unit while it is wet, it may cause trouble or you may suffer electric shock.
- ⚠ Do not poke the metal mesh with a cotton swab or a pin. The unit may break and in case it breaks it cannot be used.
- ⚠ Do not drop or apply strong shock to medicine cup. It may break and in case it breaks it cannot be used since you may suffer electric shock.
- ⚠ In order to ensure the correct use of this product, please read this manual carefully and operate in strict accordance with the requirements of the manual.
If you have any questions, please contact the supplier or manufacturer.
- ⚠ When the storage and transportation temperature is lower than 5°C, the device should be restored to normal working conditions before use.


4. Using Instructions

1. Connect the power supply:

The equipment supports external power supply (5V 1A).


Please note that this main unit is USB connector and power supply is 5V 1A. Plug the USB connector of the main unit into the power bank.

The product includes a Type-C adapter head, so you can use a power supply with a Type-C adapter head connector.


2. Fill the medication:

Open the medicine cup, fill the medication, close the medicine cup tightly. When using high concentration, high viscosity, suspension and volatile liquid, be sure to follow the doctor's instructions.

3. Connect medicine cup:

Insert the other end of the main unit connection cable into the input port of the medicine cup and make sure it is connected tightly. Do not shake violently when the medicine cup is filled with liquid.

4. Install the inhalation mask:


5. Place the mask on the face:

6. Start the inhalation:

Please make sure that the power is turned on and meets the output capacity of 5V 1A. Short press the power button of the main unit to start inhalation.

- Take a slow, deep breath and inhale the medicine.
- For agents with high surface activity or high viscosity, the nebulization rate may be reduced.
- It is normal that some sound might be produced by the metal mesh during nebulization.
- If the solution has high surface activity, a small amount of the solution may exudate the metal mesh.
- If the liquid in the medicine cup runs out, press the switch button on the main unit to turn off the device to avoid potential damage to the metal mesh.
- When the medication is about to run out, users are advised to tilt the medicine cup slightly to themselves to facilitate the remaining liquid to touch the metal mesh for nebulization.
- When inhaling with a mask, do not cover the mask vent with your hands or other objects.
- There is a vent on the lid of the medicine cup, please do not cover it with your hands or other objects, so as not to affect the normal spray.
- During the inhalation, if the user feels uncomfortable, please turn off the machine immediately and consult the doctor in time.



7. Stop the inhalation:

Short press the power button of the main unit to stop inhalation. Turn off the external power output and unplug the main unit from the external power supply.

5. Cleaning and disinfection

Clean the medicine cup and mask

- 1) Remove inhalation mask and medicine cup from the main unit.
- 2) Open the medicine cup and discard the medication.
- 3) Add a small amount (2-5ml) of water to the medicine cup, then fasten the lid tightly. Shake the medicine cup gently so that the residual liquid is fully dissolved in water. Open the medicine cup, pour the solution out of the cup.
- 4) Wash the mask in the water.
- 5) Immediately after washing the parts, please dry thoroughly in a clean environment.
- 6) Store all parts of the nebulizer in a clean place.

Disinfect the medicine cup and mask

- 1) The medicine cup can be immersed in 75% medical alcohol for more than 3 minutes, then take out the medicine cup and place it in a clean place to dry, dry thoroughly and store in a clean place.
- 2) Use 75% medical alcohol to scrub the mask for more than 3 minutes, then dry the mask and store it in a clean place.

- 3) After alcohol disinfection, rinse thoroughly with water to avoid alcohol residue. Be sure to disinfect the medicine cup and mask after each use, and replace the parts in time if the parts are seriously damaged.

Replace the medicine cup and mask

The medicine cup should be replaced every 3 months. Otherwise, safety and health requirements, functional capacity and safe treatment will no longer be guaranteed. It also needs to be replaced with a new medicine cup if the medicine cup is damaged or if the patient cannot inhale normally.

6. Maintenance

- Before each use, check the device. Do not use the device if it is damaged in any way. The continuous use of a damaged unit may cause injury, improper results, or serious danger.
- Store and use the device in a cool, dry and ventilated environment. Keep away from fire and heat sources, or the battery could explode.
- If you have any problems with this device, such as setting up, maintaining or using it, please contact your distributor or Vapo Healthcare Ca., Ltd. Don't open or repair the device by yourself.
- Please contact your distributor or manufacturer if any unexpected operation or event occurs.

7. Trouble shooting

Trouble	Possible cause	How to correct
The nebulization rate is extremely low	The metal mesh is stained	Clean the medicine cup to remove the stain. If the problem cannot be corrected by the above-mentioned procedure, replace the medicine cup with a new one.
The unit does not nebulize	The main unit is not connected correctly	Remove the main unit from the power bank and connect it correctly.
The light of main unit works correctly but the unit does not nebulize	main unit connection cable failure	Replace it with a new main unit
	The metal mesh is extremely stained	Clean the medicine cup according to the instructions. If there is no troubleshooting after cleaning, please replace the medicine cup with a new one.
Liquid leakage	The medicine cup is damaged	Replace it with a new medicine cup.

8. Disposal

DO NOT dispose of the batteries in household waste. Dispose of the batteries according to the local regulations dealing with the disposal of these special materials (e. g. to the collecting points).

The device is made of a combination of plastic and stainless steel. Be sure to dispose of it in accordance with local regulations as unsorted municipal waste. You may recycle it at your local community or appliance recycling center.

9. Electromagnetic compatibility

Guidance and Manufacturer's Declaration

- 1) This equipment needs to be installed and put into service in accordance with the information provided in the ACCOMPANYING DOCUMENTS;
- 2) IMMUNITY TEST LEVELS for basic safety and essential performance of ME equipment and ME systems should be chosen based on a high probability of maintaining basic safety and essential performance, and shall be according to the professional Healthcare facility environment, home Healthcare environment, and special environment, based on the locations of intended use.
- 3) HOME HEALTHCARE ENVIRONMENT is a dwelling place in which a patient lives or other places where patients are present, excluding professional Healthcare facility environments where operators with medical training are continually available when patients are present. Such as schools, outdoors, homes, vehicles, hotels and guesthouses.

EXAMPLE: As indicated in Table 6 of IEC 60601-1-2: 2014 for ME EQUIPMENT, a typical cell phone with a maximum output power of 2W yields $d = 3.3\text{m}$ at an IMMUNITY LEVEL of 3V/m.

A1 Electromagnetic Emissions For all EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration - electromagnetic emission		
The VP-M10 is intended for use in the electromagnetic environment specified below; The customer or the user of the P-M10 should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The VP-M10 uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The VP-M10 is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Hammonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/-flicker emissions IEC 61000-3-3	Complies	

A2 Electromagnetic Immunity
For Home Healthcare Environment EQUIPMENT and SYSTEMS

Guidance and manufacturer's declaration - electromagnetic Immunity			
The VP-M10 is intended for use in the electromagnetic environments specified below; The customer or the user of the VP-M10 should assure that it is used in such an environment.			
Immunity test	IEC60601 test level	compliance level	Electromagnetic environment -guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2kV, ±4 kV, ±8kV, ±15kV air	±8kV contact ±2kV, ±4kV, ±8kV, ±15kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2KV for power supply lines	±2KV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-8	±1kV line to line ±2kV line to ground	±1kV line to line ±2kV line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

Voltage dips, interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	0% UT; 0,5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°	Mains power quality should be that of a typical commercial or hospital environment. If the user of the VP-M10 requires continued operation during power mains interruptions, it is recommended that the VP-M10 be powered from an uninterruptible power supply.
Conducted RF IEC 61000-4-6	3Vrms 150kHz to 80MHz	3Vrms	The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:
Radiated RF IEC 610004-3	10V/m 80MHz to 2.7GHz	10V/m	$E = \frac{6}{d} \sqrt{P}$ <p>Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.</p>

NOTE: UT is the a.c. mains voltage prior to application of the test level.

10. Warranty card

MODEL		
Warranty Period	One year from purchasing date	
Purchasing shop	Name:	Telephone:
	Address:	
Client	Name:	Telephone:
	Address:	
Requiring Record		
Date	Trouble	Service Man

Guarantee regulation:

- The Guarantee is for the main unit, not applicable for the medicine cup as the shelf-life of the medicine cup varies greatly according to usage.
- The product must be accompanied by proof of purchase, either a bill of sale or other proof supporting that the device is within warranty period.
- This warranty does not apply to a product which has been damaged as a result of improper maintenance, an accident, improper voltage supply or any other form of misuse. The warranty is also void if the owner repairs or modifies the product in any way. VAPO is not liable for any incidental or consequential damages with regard to this product. The warranty also excludes any liability other than that stated above. No other warranty is given.
- Legal rights vary from country to country. Some countries do not allow the exclusion or limitation of incidental or consequential damages, so the above limitation may not apply to you.

NEBUALDO

Easy Pro

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Mesh Nebulizer

Vapo Healthcare Co., Ltd.
Southern unit of third floor,
building B, No. 99 Yudai West Rd,
High Tech district, KunShan, Suzhou,
Jiangsu Province, China

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Share Info GmbH