NEBUALDO

Hands-Free

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Portable Mesh Nebulizer

Manual de Usuario / User manual

Aldo-Unión
www.aldo-union.com

NEBUALDO

Hands-Free

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Portable Mesh Nebulizer

Manual de Usuario / User manual

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- Thank you very much for purchasing Nebualdo Hands-Free.
- In order to make sure that this product can be used correctly, please read this manual carefully before use.
- Please keep this manual in a convenient place for easy access.
- Illustrations contained in this user manual are schematic.

1. Important safety notes

- ⚠ Before use, ensure that there is no visible damage to the device or accessories. In case of any doubt, do not use the device and contact your retailer or the specified Customer Service address.
- ⚠ Do not use medical devices or medicines containing essential oils for nebulization.
- ⚠ You should always follow the instructions of your doctor regarding the type of medication to use, the dosage, and the frequency and duration of inhalation. Only use medicines prescribed or recommended by your doctor or pharmacist.
- ⚠ The use of this product for children and people with special needs must be carried out under the correct guidance and supervision.
- ⚠ This device must only be used for the intended purposes, only for nebulization. Do not use the device for any other purpose.
- ⚠ Clean and disinfect the medication cup and accessories before using or after not having used the device for a while.

- ⚠ Please stop using the device if the components are damaged or in case they fall into the water accidentally.
- ⚠ Keep the device away from your eyes when being used, as the nebulized medication could be harmful.
- ⚠ Keep packaging material out of the reach of children (risk of suffocation).
- ⚠ Do not use any additional parts that are not recommended by the manufacturer.
- ⚠ If any serious incident occurs in relation to the device, please report it to the manufacturer and the competent authority of your member state immediately.

2. Product description

1. Product name

Nebualdo Hands-Free

2. Model no.

Air Mask II

3. Operating principle and mechanism

The operating principle of the portable mesh nebulizer consists of the rapid oscillation of the high-frequency circuit, that causes the resonant oscillation of the piezoelectric ceramic transducer, which in turn makes the micro-metal mesh vibrate making the liquid pass through it and expelling the liquid quickly in the form of numerous tiny atomized particles, that pass through the mask or mouthpiece, until they reach the patient's respiratory system, in order to achieve the purpose of the inhalation therapy. The respiratory system is an open system. The liquid is atomized into particles and when the patient inhales this nebulization of the medicine, the medicine can be absorbed directly once

deposited in the patient's mouth, throat, trachea, bronchi, alveoli, etc., as it is absorbed through the tissue of the mucous membrane in order to achieve the purpose of the treatment.

4. Applicable scope and intended user


The device is a mesh nebulizer, designed to aerosolize liquid medicines for inhalation therapy in professional healthcare environment and in home healthcare environment. Suitable for pediatric and adult patient, infants, children and people with disabilities. It should be used under supervision of an adult.

5. Specification

Power Supply	DC 3.7V (Li-ion Battery) or DC 5V 1A with IEC 60601-1 approved AC adapter
Nebulization Rate	0.15ml/min- 0.90ml/min
Particle Size	MMAD<5 μ m
Power Consumption	<4.0W
Working Frequency	130 kHz \pm 10 kHz
Medication Cup Capacity	6ml
Noise	\leq 50dB
Product Dimensions/ Weight	72.5 mm (Length) x 67.5 mm (Width) x 50.5 mm (Height)/72 g

Working Environment	Temperature, 5°C to 40°C Relative Humidity \leq 80% R.H. Non-condensing state Atmospheric pressure: (86.0- 106.0) kPa
Storage/Delivery Environment	Temperature: -20°C to 55°C Relative Humidity \leq 80% R.H. Non-condensing state Atmospheric pressure: (86.0- 106.0) kPa

The mean particle size in this nebulizer is measured with 0.9% physiological saline solution under conditions of a temperature of 25°C and a relative humidity of 59%. And the equivalent particle size distribution curve of the mist particles measured under these conditions is as follows:


Note: the horizontal axis is the particle size value, the value is logarithmic distribution; the left vertical axis is the cumulative percentage of the volume, corresponding to the rising trend of the curve; the right vertical axis is the percentage of the volume of a certain section, corresponding to the histogram's undulating graph.

6. Product composition

Nebualdo Hands-Free consists mainly of the main unit and medication cup, spray assembly and micro USB cable.

7. Product contents



A. Medication cover
B. Medication cup
C. Spray nozzle
D. Mesh diaphragm
E. Main unit
F. Power switch

G. Electrode contact pin
H. USB charging port
J. Adult mask (PVC)
I. Child mask (PVC)
K. Mouthpiece (PP)
L. Headband

Schematic diagram of host structure

3. Installation instructions

1. Remove all packings, then remove the unit and accessories .
2. Assemble the medication cup to the main unit. When you install it, you should hear the clear clasp sound (as shown in the schematic diagram of the installation of the medication cup).
3. Install the suction mask and the mouthpiece as shown in the schematic diagram.


Note: Clean, disinfect and dry all parts before installation.

4. Power Supply

- 1) The nebulizer has a USB cable for charging, since power adapter is not included, please use an IEC 60601-1 approved AC adapter (output: DC 5.0V 1.0A) for charging .
- 2) The power system of the device is equipped with one rechargeable lithium battery.
- 3) When the device is running low, please charge it with the help of the Micro USB cable, then it will able to work again .

Charging diagram

⚠ Note: Before charging, please make sure the power adapter where you are going to connect it has an DC 5.0V 1.0A output.

⚠ Note: This instrument is independent charging, please do not charge it with any other electronic equipment.

4) Battery Charging

- a. The battery group can supply power up to 60 minutes continually after full charging.
- b. When the low battery level is detected, the blue indicator light flashes for 5 times and turns it off.
- c. Please use the power adapter (DC5.0V, 1A) to charge the batteries for about 2 hours.
- d. The blue indicator light flashes while charging and keeps on when it is fully charged.

⚠ Notes:

- 1) Battery has been loaded, do not disassemble it yourself.
- 2) Rechargeable batteries must not be replaced by user, only by the manufacturer.
- 3) Charge the device at least once a month for storage periods exceeding one month.
- 4) Alkaline, lithium-ion or other batteries are not indicated for this device.
- 5) Please charge at least 30 minutes before using the device for the first time.
- 6) In order to achieve as long a battery service life as possible, fully charge the battery at least once a month .

⚠ Warning:

- Please dispose of the used batteries according to the local environmental regulation. Do not dispose of them with the household waste to avoid environment pollution.
- Do not dismantle or repair the equipment or its components, do not dismantle or replace the battery. If you need to replace the battery, please consult or purchase it from Feellife Health.
- If your skin or eyes come into contact with fluid from a rechargeable battery cell, wash the affected areas with water and seek medical assistance .
- Choking hazard! Young children may swallow and choke on batteries. Store batteries out of the reach of young children.
- Risk of explosion! Never throw batteries into a fire.
- Do not disassemble, split or crush the batteries .
- Only use chargers as specified in the instructions for use.
- Batteries must be charged correctly prior to use. The instructions from the manufacturer and the specifications regarding correct charging given in these instructions must be observed at all times.

4. Instructions for use


1. Indicator description form

Green light keeps on (uniform nebulization mode)	Working
Green light flashes (micro-wave mode)	
Green light flashes 3 times	Long press for 3 seconds to switch the atomization mode
Blue light flashes 5 times	Battery low, switch off
Orange light flashes 10 times	Without medicine and switch off
Blue light flashes	Charging
Blue light keeps on	Fully charged
Green light flashes 10 times	Set time at 10 minutes

Note: There is enough liquid, but the orange light flashes ten times, while the orange light is flashing, please press the power button once so that it continues working.

- Preparation: clean and disinfect the components of the spray mouthpiece and medication cup before use.
- Injection of the liquid: open the cap rotating it counterclockwise, inject the solution and turn the cap clockwise to close.


Note: Please add the liquid before turning on the device.


Note: Leakage protection! When pouring medication into the medication cup, make sure that you only fill it up to the maximum mark (6ml). The recommended fill volume is between 2 and 6ml. Nebulisation only occurs while the substance to be nebulised is in contact with the mesh. If this is not the case, nebulisation stops automatically. So please try to hold the device as vertically as possible.

4. Nebulization

- Connect the mask or mouthpiece, press  power button to turn it on, start the nebulization.
- It starts with a medium spray rate, and the spray rate can be adjusted by pressing “+” & “-” as needed ;
- This device has a Micro-wave mode in which the spray rate is changed automatically from max to min rate or conversely in 6 seconds. After switching to the micro-wave mode, long press  power button for 3 seconds (when the green light flashes 3 times) to switch to the uniform nebulization mode.
- For the start before the atomization, shake the level slightly, so that the solution comes into complete contact with the atomization using the three following inhalation ways according to individual needs.

- 5) Slowly take a deep breath and then inhale the drug.
- 6) The nebulizer will switch off automatically after 10 minutes. If you need to continue using it, press the Power Button . Please make sure there is enough liquid in the medication cup.
- 7) After nebulization, press the Power Button  to turn off the nebulization. Pour out the trace of residual liquid into the medication cup and do not reuse (disassemble as shown).
- 8) Use purified water to clean the medication cup, spray assembly, lid, thread set and accessories, and then follow the recommended method of disinfection.

Notes:

- 1) Hold on the nebulizer to keep the liquid and nebulization mesh disc in full contact, a slight swing does not affect the use but do not lean back the nebulization.
- 2) Receive the treatment according to the doctors' recommendation. Keep quiet and relaxed in the course of the treatment.
- 3) In case the liquid gets coagulated around the spray assembly and mesh disc, which affects the nebulization effect after nebulization. You must stop the nebulization and remove the mouthpiece and other accessories. Then use a clean medical gauze to wipe the residue. Do not touch the mesh disc centre spray area to prevent damage to the mesh disc.
- 4) Essential oils, cough syrups, gargling solutions and drops to be used as a rub or in a steam bath are totally unsuitable for inhalation using a nebulizer. These additives are often viscous and can impair the correct functioning of the device and therefore, the effectiveness of the application in the long term.

5. Cleaning and disinfection method

After each use, it is necessary to clean and disinfect the cup components (including medication cup and the lid), spray assembly, mask or mouthpiece. The recommended specific cleaning and disinfection methods are the following:

1. Cleaning

Please turn off the power when cleaning machine. Do not connect the machine to power supply.

1) Remove the components from the main unit: the medicine cup components (including medication cup and the lid), spray assembly, mask or mouthpiece; and soak all the components in clean warm water (below 40 degrees Celsius) for about 5 minutes.

2) After cleaning, wipe all the components with a clean and sterile medical gauze, and keep them sufficiently dry.

3) Wipe the outer shell of the main unit. If there is residue medicine remaining at the electrode contact, please clean it with a wet sterile medical gauze. After cleaning, keep the main unit dry.

4) Store all parts in a dry and clean place to avoid contamination.

Notes:

- The main unit must not be washed with water to prevent water from entering the main unit.
- Use a clean sterile gauze to wipe off water on the main unit and components and keep them dry to ensure a safe use next time.
- The masks must not be placed in hot water!

2. Disinfection

After each use, it is necessary to disinfect the components of the cup, including medication cup, lid, spray assembly, mask or mouthpiece, etc., as follows:

1) Disinfection with hydrogen peroxide

Disinfect all the components by placing them in 3% hydrogen peroxide for 10 minutes, including cup components (including medication cup and the lid), spray assembly, mask or mouthpiece. Or wipe and disin-

fect the medicine cup component (including medication cup and the lid), spray assembly, mask or mouthpiece. After disinfection, rinse all the parts with distilled water, then wipe with a clean and sterile medical gauze or air-dry naturally to keep all parts dry.

A. Please study the instructions for use of hydrogen peroxide and do not immerse in this solution for a long time.

B. Do not use strong oxidizing agents such as perchlorate or disinfectants that are corrosive to metals, polymer compounds or polymers.

2) Disinfection with ethanol

Place the liquid cup components (including medication cup and the lid), spray assembly and mask or mouthpiece in 75% medical ethanol for 10 minutes for disinfection. Or wipe and disinfect the medicine cup component (including medication cup and the lid), spray assembly, mask or mouthpiece. After disinfection, rinse all the parts with distilled water, then wipe with clean and sterile medical gauze or air-dry naturally to keep all parts dry.

Note:

- Disinfectants remaining on components need to be wiped off with a sterile medical gauze to ensure a safe use next time. Do not touch the central area of the mesh diaphragm when washing or cleaning the spray connector, so as to avoid damaging the mesh diaphragm.

3) Drying

- Dry the parts carefully using a soft cloth.
- Shake the spray assembly gently from side to side (5-10 times), so that the water inside the mesh is removed from the tiny holes.

- Place the individual parts on a clean, dry and absorbent surface and allow to dry completely (at least 4 hours). Note: Please ensure that the parts are completely dry after cleaning, otherwise the risk of bacterial growth increases. Put the parts together again if they are completely dry and place the parts in a dry, sealed container. Make sure that the spray assembly is completely dried off with the shaking. Otherwise, the nebulization may not work after reassembling the device. If this is the case, shake the spray assembly again so that the water can escape. The nebulizer should then work as normal.

6. Storage and maintenance

1. Nebulizer storage

1) Storage conditions:

- a. Room temperature: -20°C to 55°C
- b. Relative humidity: $\leq 80\%$ R.H. Non-condensing state
- c. Atmospheric pressure: 86.0 kPa-106.0 kPa
- d. Others: Non-corrosive gas, good ventilation, avoid high temperature, humidity and direct sunlight.

2) Storage instructions:

- a. The device is valid for 5 years in the above-mentioned storage conditions.
- b. The nebulizer should be promptly cleaned and disinfected after use, the medication cup and other accessories should be stored into the packing box once they are completely dry. They should be stored in the required environmental condition, avoiding impacts.

2. Nebulizer maintenance

1) Normal working conditions

- a. Room temperature: 5 °C to 40 °C
- b. Relative humidity: $\leq 80\%$ R.H. Non-condensing state
- c. Atmospheric pressure: 86.0-106.0 kPa
- d. Power source: DC 3.7V (Li-ion battery integrated) or DC 5.0V with approved AC adapter (purchased by user, for technical specification please refer to Chapter 2 for further details).

2) Maintenance instructions:

- a. Please use nebulizer under normal conditions;
- b. Do not use the nebulizer near the heating device or open flame. Do not use a microwave oven, fan or other devices in order to dry nebulizers and accessories.
- c. Do not expose the nebulizer and accessories to corrosive liquids and gases.
- d. Do not wrap the power cable around the unit.
- e. When using a nebulizer, if any irregularities are encountered, seek a solution in accordance with Chapter 8.
- f. Dry the parts immediately after washing. Store the device and the components in an environment that meets the requirements, be careful to avoid collisions.
- g. Direct sunlight, lint and dust may cause vibrating mesh getting rusted and oxidized and decrease nebulization rate.
- h. If the nebulizer still does not work properly, please contact the manufacturer or the distributor.

7. Contraindications, precautions, notices and warnings

1. Contraindications

- 1) This product is not suitable for Pentamidine drugs.
- 2) The use of this product is prohibited in patients with pulmonary oedema.
- 3) The use of this product is prohibited in patients with acute asthma and acute pulmonary infarction episodes.

2. Notices and warnings

<ul style="list-style-type: none"> • The nebulizer is a medical device and it is intended for human use only. Please follow the instructions in the manual or use it under the guidance of a doctor. Infants and children and people with special needs should use it under the supervision of their guardian. 	
<ul style="list-style-type: none"> • Please use original parts and accessories. Warranty service is not provided for damage caused by accessories not included in our list and damaged due to user's personal causes. 	
<ul style="list-style-type: none"> • The equipment waterproof rating is IP22, the main unit cannot be washed, and the ingress of water should be prevented. 	
<ul style="list-style-type: none"> • Please refer to the "Troubleshooting" section and contact the service centre for maintenance. Do not attempt to repair the equipment yourself. 	
<ul style="list-style-type: none"> • Please clean and disinfect the unit after use and refer to the clean and disinfect section. 	
<ul style="list-style-type: none"> • The nebulizer is for medication atomization 	

<ul style="list-style-type: none"> • Remove the liquid before storing the device. 	
<ul style="list-style-type: none"> • Please verify that all the accessories are intact before use. 	
<ul style="list-style-type: none"> • The mouthpiece is a disposable accessory, only for person single-use to avoid cross-infection. Max reusable time up to 3 times. 	
<ul style="list-style-type: none"> • Try to make the liquid fully reach the mesh disc when you use the device. 	
<ul style="list-style-type: none"> • Never submerge the device in water and do not use it in the bathroom. Under no circumstances should liquid enter the device. 	
<ul style="list-style-type: none"> • Do not use the unit near flammable gas atmospheres or near oxygen and anaesthetic mixtures. 	
<ul style="list-style-type: none"> • The device is expected to be used in an electromagnetic environment with controlled radioactive radiation, as far away as possible from the radiation source. 	
<ul style="list-style-type: none"> • Do not use the machine in a high-temperature environment as it may cause problems and fire. 	
<ul style="list-style-type: none"> • Keep the device and parts away from strong vibrations, such as impacts. 	
<ul style="list-style-type: none"> • Please do not use liquids which contain esters, oil or suspended particles, including herbal extract. 	
<ul style="list-style-type: none"> • Do not wash the whole unit under running water to avoid water come inside the device, especially inside the USB connector. 	
<ul style="list-style-type: none"> • Do not use the microwave ovens to dry or disinfect the unit as it could cause a fire. 	
<ul style="list-style-type: none"> • Keep it out of the reach of children, infants and psychotic patients. 	

<ul style="list-style-type: none"> • Please use lithium battery power chargers (output 5.0V/1.0A) from qualified manufacturers. 	
<ul style="list-style-type: none"> • Do not touch the centre of spray mesh with your hand or with other sharp objects. It may cause damage, and it could not be used. 	
<ul style="list-style-type: none"> • The nebulizer has no moisture-proof or dust-proof function, and the product should not be stored in a wet or dusty environment. 	
<ul style="list-style-type: none"> • If the device has been dropped, exposed to high levels of moisture or suffered any other damage, it must no longer be used. In case of doubt, contact Customer Service or the distributor. 	
<ul style="list-style-type: none"> • Do not disassemble, repair or modify the device as it may cause electric shock, leakage or fire. 	
<ul style="list-style-type: none"> • Dispose of the wasted device and accessories in accordance with local regulations. Illegal disposal may cause pollute environment. Users can dispose of the device in an appropriate manner in order to avoid environment pollution. 	

3. Warnings

<ul style="list-style-type: none"> • Please consult your doctor before using the device if you have diabetes or other illnesses. 	
<ul style="list-style-type: none"> • The use and purchase of the device should be advised by a doctor, please refer to the doctor's advice regarding the medication type, dosing and way of use. 	

<ul style="list-style-type: none"> • Please stop using the device if you feel uncomfortable and seek medical help. 	
<ul style="list-style-type: none"> • Volatile oil is not allowed, may cause damage to the module. 	
<ul style="list-style-type: none"> • Only water-soluble medicines containing alcohol or saline-diluted medicines can be used for nebulization treatment, otherwise bronchospasm may be caused. 	
<ul style="list-style-type: none"> • Oily medication is not allowed. 	
<ul style="list-style-type: none"> • The device is not compatible with inhalational anaesthesia systems and breathing systems. 	
<ul style="list-style-type: none"> • Check the product leaflet for any contraindications for use with the usual systems for aerosol therapy. 	
<ul style="list-style-type: none"> • Do not use any liquids with a viscosity greater than 5, as this can cause irreparable damage to the mesh. 	
<ul style="list-style-type: none"> • Only use the medication with a saline solution. 	
<ul style="list-style-type: none"> • When the user needs to carry out product maintenance, the manufacturer can provide the product circuit diagram, key component list and other information according to the requirements. For further details or advice on maintenance, the user can contact the qualified technical staff of the after-sales service of the manufacturer. 	
<ul style="list-style-type: none"> • Do not spill liquid on the device to avoid leakage, the possibility of electric shock and causing malfunction or failure to use. 	

8. Troubleshooting TIPS

Troubleshooting for nebulizer		
Item	Trouble	Possible cause / solution
1	It doesn't work when turned on	1) Make sure the nebulizer has enough power, when the blue and green indicator light are always on, please recharge the device. 2) Check if the medication cup is fully filled.
2	Low volume	1) Check if the medication cup has been filled with the right medication, which should be water-soluble and non-corrosive. 2) Check if the amount of liquid is too high. 3) The spray mesh could be blocked, you can drop 2 or 3 drops of white vinegar into the medication cup with 3-6 ml of water, then fully nebulize it. Clean the medication cup to disinfect it for next nebulization.
3	Which medication is most suitable for nebulization	1) Please consult a doctor and follow doctor's advice. 2) Don't use with sticky medication.
4	Some medication residual	It is a normal sign, if there are some strange noises or if the unit turns off due to insufficient medication, please stop nebulization.

5	Special care for babies and children	1) For infants, the mask should cover the nose and mouth to ensure effective inhalation. 2) For children, the mask should also cover both the nose and mouth. It is not very useful to use the nebulizer on someone who is sleeping, as not enough medication will reach the lungs in this case. Note: Inhalation should only be carried out under the supervision of an adult and with their assistance, and the child should not be left alone.
6	Each user needs an individual consumable item	Each user should use individual consumable items, including mask and mouthpiece.

9. Disposal

Battery disposal

The empty, completely flat rechargeable batteries must be disposed of through specially designated collection bins, recycling points or electronics retailers. You are legally required to dispose of the batteries.

General disposal

For environmental reasons, do not dispose of the device with the household waste at the end of its useful life. Dispose of the device at a suitable local collection or recycling point in your country. Dispose of the device in accordance with EC Directive - WEEE (Waste Electrical and

Electronic Equipment). If you have any questions, please contact the local authorities responsible for waste disposal.

10. After-sales service

1. Please contact our after-sales service department to obtain warranty service.
2. If necessary, you can provide the circuit diagram and the necessary information for the repair. If there is any problem in the maintenance of the electrical circuit, you can contact the manufacturer.


For further information please refer to the user manual.

*The company reserves the right of final interpretation of the warranty card, which may be subject to change without prior notice.

11. Symbol description & electromagnetic compatibility

1. Signs and symbols

 <p>Refer to instruction manual/ booklet NOTE ON ME EQUIPMENT</p>	 <p>Separate collection for electrical and electronic equipment</p>	 <p>CE identification+ notified body number</p>
--	--	--

 <p>General imperative</p>	 <p>ON/OFF Switch</p>	 <p>keep away from sunlight</p>
<p>IP22</p> <p>Waterproof Grade</p>	 <p>Fragile</p>	 <p>keep dry</p>
 <p>The application part of Type BF</p>	 <p>Note, Warning, refer to enclosed file</p>	 <p>Cause for electric shock</p>
 <p>Decrease nebulization rate</p>	 <p>Increase nebulization rate</p>	 <p>Manufacturer</p>
 <p>Batch code</p>	 <p>Manufacturing date</p>	 <p>Charging Interface</p>
 <p>Precautions</p>	 <p>Authorized representative in the European Community</p>	 <p>Medical Device</p>

2. EMC declarations

<ul style="list-style-type: none">• Nebualdo Hands-Free Portable Mesh Nebulizer meets the requirements of electromagnetic compatibility in IEC60601-1-2.
<ul style="list-style-type: none">• The user must install it and use it according to the attached electromagnetism compatibility information.
<ul style="list-style-type: none">• Portable and mobile RF communication devices and some household appliances, such as mobile phone, interphone, microwave oven, dry blower, may influence, so the nebulizer should be kept away from them when in use.
<ul style="list-style-type: none">• Guidance and manufacture's declaration stated in the appendix.

3. EMC information.

With the increased number of electronic devices such as PC's and mobile (cellular) telephones, medical devices in use may be susceptible to electromagnetic interference from other devices. Electromagnetic interference may result in incorrect operation of the medical device and create a potentially unsafe situation.

Medical devices should also not interfere with other devices.

In order to regulate the requirements for EMC (Electro Magnetic Compatibility) with the aim of preventing unsafe product situations, the IEC60601-1-2 standard has been implemented. This standard defines the levels of immunity to electromagnetic interferences as well as maximum levels of electromagnetic emissions for medical devices.

Medical devices manufactured by FEELLIFE HEALTH INC. conform to this IEC60601-1-2 standard for both immunity and emissions.

Nevertheless, special precautions need to be observed:

1. NOTE: The EMISSION characteristics of this equipment make it suitable for use in a residential environment (for which CISPR 11 class B is normally required)

2. WARNING: The use of accessories and cables other than those specified by FEELLIFE, with the exception of cables sold by FEELLIFE as replacement parts for internal components, may result in increased emission or decreased immunity of the device.

3. WARNING: Medical devices should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the medical device should be observed to verify normal operation in the configuration in which it will be used.

4. MOBILE RF communication equipment (including peripherals such as antenna cables and external antennas) should not be used closer than 30cm (12 inches) to any part of [ME EQUIPMENT or ME SYSTEM], including cables specified by the MANUFACTURER. Otherwise, degradation of the performance of this equipment could occur.

5. Do not use mobile telephones and other devices (such as MRI, diathermy, electrocautery, RFID and electromagnetic security systems) which generate strong electrical or electromagnetic fields, near the medical device. This may result in an incorrect operation of the unit and cause a potentially unsafe situation. Recommendation is to keep a minimum distance of 7m. Verify that the current operation of the device is correct in case the distance is shorter.

Cables	Length (m)	Signal shielding	Note
Cable USB	1	No	-

Guidance and Manufacturer's declaration - electromagnetic emissions		
This device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.		
Electromagnetic emission IEC 60601-1-2		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions C1SPR11	Group 1	This device uses RF energy only for its internal function. Therefore, its RF emissions are very low and they are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	This device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's declaration - electromagnetic immunity			
The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment- guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8kV contact ±2 kV,±4 kV, ±8 kV, ±15 kV, air	±8kV contact ±2 kV,±4 kV, ±8 kV, ±15 kV, air	Floors should be made of wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Main power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±0,5 ±1kV line to line ±0,5 kV ±1 kV, ±2 kV line to earth	±0,5 ±1 kV, line to line	Main power quality should be that of a typical commercial or hospital environment.

1. Field strengths from fixed transmitters, such as base stations for radio, telephones (cellular/cordless) and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the above-mentioned applicable RF compliance level, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating this device.

2. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the device

The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help to prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.

Maximum output power rate of transmitters (w)	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz d=1.2√P	80 MHz to 800 MHz d=1.2√P	800 MHz to 2,5 GHz d=2.3√P
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

⚠ PRECAUTIONS

1) In order to comply with EMC with requirements the aim of preventing unsafe product situations the EN60601-1-2 standard has been implemented. The nebulizer Nebualdo Hands-Free conforms to the EN60601-1-2:2015 standard for both immunity and emissions.

2) Portable and mobile RF communication equipment may affect the performance of the nebulizer, avoid strong electromagnetic interference when using it such as getting close to mobile phones, microwave ovens and so on.

⚠ WARNING

1) Equipment or systems should not be used adjacent to or stacked with other

equipment, if they must be used close to or stacked, it should be observed that the verification in its use of the configuration can be normal operation.

2) When the manufacturer of the unit sells the internal components, accessories and cables as spare parts, it is expected an increase in the emissions of this nebulizer or a reduction in immunity.

12. Configuration list

Item	Quantity	Included	
		YES	NO
Mesh Nebulizer device	1	√	
Adult Mask	1	√	
Child Mask	1	√	
Additional Mesh	1	√	
Mouthpiece	1	√	
Micro-USB Charging Cable	1	√	
Leaflet	1	√	
User Manual	1	√	
Hard Case	1	√	

Warranty Card

Check for distributor

STAMP

Type of product: _____

Serial No: _____

Name of costumer: _____

Date of purchase: _____

Contact number: _____

Address: _____

Sales Company: _____

Contact number: _____

Sales Company Address: _____

E-mail: doris_liu@feellife.com

**This warranty card is made in duplicate.
This copy is to be kept by distributor.**

- The period of free service for maintenance service for the main unit is 1 year. For the medication cup is 6 months.
 - During the warranty period, we will decide whether to repair or replace damaged parts and accessories according to the case.
 - When the warranty period is exceeded, the corresponding charge must be made according to the case.
 - The accessories (mask, mouthpiece) do not come in the warranty category since they are consumable items. For details, please refer to the manual.
 - The following cases are not included in the warranty.
 - When damage has been caused by human behavior like holding, touching, soaking and wetting.
 - When damage has been caused by not operating the product correctly according to the manual.
 - When damage has been caused by accident.
 - When the user disassemble and change the product without authorization.
 - When there is no invoice, warranty card, the product serial is broken and cannot be identified.
- Attention:** please contact the local dealer or our company with the customers tub "of the warranty card" when the machine needs maintenance; please keep the packaging of product so that it could be used when the machine is returned to maintenance.

Check for the user

STAMP

Type of product: _____

Serial No: _____

Name of costumer: _____

Date of purchase: _____

Contact number: _____

Address: _____

Sales Company: _____

Contact number: _____

Sales Company Address: _____

E-mail: doris_liu@feellife.com

**This warranty card is made in duplicate.
This copy is to be kept by the user and shown at warranty's time.**

- The period of free service for maintenance service for the main unit is 1 year. For the medication cup is 6 months.
 - During the warranty period, we will decide whether to repair or replace damaged parts and accessories according to the case.
 - When the warranty period is exceeded, the corresponding charge must be made according to the case.
 - The accessories (mask, mouthpiece) do not come in the warranty category since they are consumable items. For details, please refer to the manual.
 - The following cases are not included in the warranty.
 - When damage has been caused by human behavior like holding, touching, soaking and wetting.
 - When damage has been caused by not operating the product correctly according to the manual.
 - When damage has been caused by accident.
 - When the user disassemble and change the product without authorization.
 - When there is no invoice, warranty card, the product serial is broken and cannot be identified.
- Attention:** please contact the local dealer or our company with the customers tub "of the warranty card" when the machine needs maintenance; please keep the packaging of product so that it could be used when the machine is returned to maintenance.

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Lotus NL B.V.Koningin Julianaplein 10, 1e Verd,